



Bluetooth  
**CARDIOLINE**

## ECG



Configurator



**Mortara**



# ECG



## Important note

The Cardioline HD+ ECG can be used as a stand alone PC-ECG unit.  
The Cardiosoft option only includes the Expair software interface with Cardiosoft but no ECG module.  
The Mortara option includes the Xscribe CPI software and the Mortara ECG module.  
The Custo-Med option only includes the Expair software interface with Custo-Med but no ECG module.

Start with the base module and follow this path :

- 1 ECG
- 2 Computer
- 3 Trolley

## 1 ECG

- ECG wireless 12 leads Cardioline HD+ (rest and stress) Ref: ECG-00-01
  - Suction electrodes DT 100 (mobile version) Ref: ECG-00-02
  - Non invasive Blood pressure (NIBP) Tango external module Ref: BP-00-04
- Mortara ECG
  - ECG wired 12 leads Mortara Xscribe with patient cable IEC type (european) Ref: ECG-00-06
    - Non invasive Blood pressure (NIBP) Tango external module Ref: BP-00-03
  - ECG wireless 12 leads Mortara Xscribe with patient cable IEC type (european) Ref: ECG-00-07
    - Non invasive Blood pressure (NIBP) Tango external module Ref: BP-00-03
  - ECG wired 12 leads Mortara Xscribe with patient cable AHA type (american) Ref: ECG-00-10
    - Non invasive Blood pressure (NIBP) Tango external module Ref: BP-00-03
  - ECG wireless 12 leads Mortara Xscribe with patient cable AHA type (american) Ref: ECG-00-11
    - Non invasive Blood pressure (NIBP) Tango external module Ref: BP-00-03
- Expair / ECG GE Cardiosoft interface Ref: ECG-00-05
  - Non invasive Blood pressure (NIBP) Tango external module Ref: BP-00-02
- Expair / ECG Custo-Med interface Ref: ECG-00-14

## 2 Computer See IT Hardware configurator

## 3 Trolley See Trolley configurator

All the options listed are available depending of the products and may be unavailable within some countries.